# kementec

# TMB-D Blotting

Precipitating HRP Substrate Cat. No. 4600

### **Product Characteristics**

TMB (3, 3', 5, 5'-tetramethylbenzidine) is a chromogenic substrate for Horseradish Peroxidase (HRP). TMB produces a blue color during the enzymatic degradation of hydrogen peroxide by HRP.

TMB-D Blotting (liquid, ready-to-use substrate) is a highly active and stable blotting substrate utilized for measuring HRP probe activity. A stable blue precipitate is formed at the reaction site.

The substrate does not contain NMP (1-methyl-2pyrrolidone) making it REACH Restricted Substances List Annex XVII compliant, while ensuring maximal safety during use, and minimal waste problems after use.

#### **Composition & Properties**

Ready-to-use substrate: Includes substrate buffer and hydrogen peroxide. No other reagents should be added.

## **Working Procedure**

The following procedure is applicable to nitrocellulose membranes. The procedure must be optimized for other membranes.

- The desired amount of substrate is poured into a sealed container and allowed to reach room temperature, in the dark, before use.
- 2. After the last incubation with HRP-labelled Streptavidin or HRP-labelled secondary antibody it is recommended to wash the membrane in a 0.1 M Tris buffer pH 7.4.

- 3. Shake off the excess buffer and incubate the membrane in the TMB-D Blotting solution for 10 minutes.
- 4. Wash the membrane in distilled water and allow it to dry.

CHROMOGENIC SUBSTRATES

5. The site of positive reaction will appear blue with no or very little background staining.

## **Tips & Tricks**

- The membrane can be blocked with Kementec's Synthetic Blocking Buffer for Blotting, (cat. no. 4650).
- For long-term preservation of the results, the membranes must be stored in the dark.

## Handling & Storage

- Store solution at 2-8°C in the dark.
- Avoid exposure to light, heat and contamination with metal ions or peroxidase.
- Re-dispense only into bottles made of High-Density Polyethylene (HDPE), amber color.
  Dispensing guidelines are available upon request.

Our Immunoassay Solutions are eco-friendly, creating a healthy work environment and preserving natural resources, while helping our customers fulfill significant regulatory requirements. **We call it ECO-TEK.** 

Kuldyssen 10 DK-2630 Taastrup Denmark VAT: DK 25 80 67 35

| Tel: +45 39 27 17 71 | Fax: +45 39 20 01 78 support@kementec.com www.kementec.com