

Antibody Enhancer

BUFFERS & STABILIZERS

AP/HRP Conjugate Stabilizer/Amplifier

Cat. No. 4660

Lot No.: xxxxxxXxx**Date of manufacture:** xx-xxxx**Expiry date:** xx-xxxx (When stored at 2-8 °C in an unopened bottle)**Analysis:** The product is tested by measuring the pH in the final solution and by visual inspection of the color.**Specifications:** pH: 5.50 – 5.90
Color: Clear solution with no distinct color.**Results:** pH: x.xx
Color: xxxxxx.....**Storage:** Short-term storage: Ambient.
Long-term storage: 2-8 °C.

It is hereby confirmed that the product was manufactured according to our standard operating procedures and that the product complies with all specifications.

Performed by: Date: xx-XXX-xxxx Signature: XX**Approved by:** Date: xx-XXX-xxxx Signature: XX

Our Immunoassay Solutions are eco-friendly, creating a healthy work environment and preserving natural resources, while helping our customers fulfill significant regulatory requirements. **We call it ECO-TEK.**