

# Protein-StabilPLUS

## BUFFERS &amp; STABILIZERS

**Antibody/Antigen Stabilizer**

Cat. no. 4720

**Product Characteristics**

Protein-StabilPLUS is used for the dilution and stabilization of antibodies and other proteins and biomolecules in solution. This enables the making of pre-diluted, ready-to-use standards and controls, minimizing assay errors generated when dissolving lyophilized standards.

**Composition & Properties**

Protein-StabilPLUS is a ready-to use-buffer that appears clear. The product is based on a mild Tris buffer containing proprietary stabilizing components. Protein-StabilPLUS contains neither BSA, nor other material from bovine serum, azide, mercury or other toxic components.

**Working Procedure**

1. Make a series of dilutions of the protein in Protein-StabilPLUS in order to determine the optimal dilution factors.
2. Run the assay as usual or store the diluted protein preferably at 2-8 °C

**Handling & Storage**

- Store solution at 2-8 °C.



Our Immunoassay Solutions are eco-friendly, creating a healthy work environment and preserving natural resources, while helping our customers fulfill significant regulatory requirements. **We call it ECO-TEK.**