

# Effect Diluent High

## BUFFERS & STABILIZERS

**Sample Diluent & RF/HAMA Blocker**

Cat. No. 7090

**Lot No.:** xxxxxxXxx

**Date of manufacture:** mm-yyyy

**Expiry date:** mm-yyyy (When stored at 2-8 °C in an unopened bottle)

**Analysis:** The diluent is inspected visually for appearance and the pH is measured.

**Specifications:** Appearance: Solution with slightly yellow to no distinct color.  
pH: 7.00 – 7.40

**Results:** Appearance: Passed  
pH: Passed

**Storage:** Short-term storage: Ambient.  
Long-term storage: 2-8 °C.

**Performed by:** **Date:** dd-MMM-yyyy **Signature:** XX

**Approved by:** **Date:** dd-MMM-yyyy **Signature:** XX

Our Immunoassay Solutions are developed in order to create a healthy and safe work environment and preserving natural resources, while helping our customers fulfill significant regulatory requirements.  
**We call it ECO-TEK!**



DOC ID: 04-5-05-422 Ver. 4