

## UNI-StabilPLUS

## BUFFERS &amp; STABILIZERS

## AP/HRP Conjugate Stabilizer

Cat. No. 5230

Lot No.: xxxxxxXxxDate of manufacture: xx-xxxxExpiry date: xx-xxxx (When stored at 2-8 °C in an unopened bottle)

**Analysis:** The product is tested by measuring the pH in the final solution and by visual inspection of the color.

**Specifications:** pH: 6.60 – 7.00  
Color: Clear solution with no distinct color.

**Results:** pH: x.xx  
Color: xxxxxx.....

**Storage:** Short-term storage: Ambient.  
Long-term storage: 2-8 °C.

It is hereby confirmed that the product was manufactured according to our standard operating procedures and that the product complies with all specifications.

**Performed by:** Date: xx-XXX-xxxx Signature: XX

**Approved by:** Date: xx-XXX-xxxx Signature: XX



Our Immunoassay Solutions are eco-friendly, creating a healthy work environment and preserving natural resources, while helping our customers fulfill significant regulatory requirements. **We call it ECO-TEK.**