

Sample X-tra

BUFFERS & STABILIZERS

Sample Dilution Buffer

Cat. No. 5350

Lot No.: xxxxxxXxxDate of manufacture: xx-xxxxExpiry date: xx-xxxx (When stored at 2-8 °C in an unopened bottle)

Analysis: The product is tested by measuring the pH in the final solution and by visual inspection of the color.

Specifications: pH: 6.70 – 7.10
Color: Clear solution with a slightly yellow to no distinct color

Results: pH: x.xx
Color: xxxxxx.....

Storage: Short-term storage: Ambient.
Long-term storage: 2-8 °C.

It is hereby confirmed that the product was manufactured according to our standard operating procedures and that the product complies with all specifications.

Performed by: Date: xx-XXX-xxxx Signature: XX

Approved by: Date: xx-XXX-xxxx Signature: XX



Our Immunoassay Solutions are eco-friendly, creating a healthy work environment and preserving natural resources, while helping our customers fulfill significant regulatory requirements. **We call it ECO-TEK.**