

Effect Diluent Medium

BUFFERS & STABILIZERS

Sample Diluent & RF/HAMA Blocker

Cat. No. 7080

Lot No.: xxxxxxXxx

Date of manufacture: mm-yyyy

Expiry date: mm-yyyy (When stored at 2-8 °C in an unopened bottle)

Analysis: The diluent is inspected visually for appearance and the pH is measured.

Specifications: Appearance: Solution with slightly yellow to no distinct color.
pH: 7.00 – 7.40

Results: Appearance: Passed
pH: Passed

Storage: Short-term storage: Ambient.
Long-term storage: 2-8 °C.

Performed by: **Date:** dd-MMM-yyyy **Signature:** XX

Approved by: **Date:** dd-MMM-yyyy **Signature:** XX

Our Immunoassay Solutions are developed in order to create a healthy and safe work environment and preserving natural resources, while helping our customers fulfill significant regulatory requirements.
We call it ECO-TEK!



DOC ID: 04-5-05-421 Ver. 4